

BATTERY INSTALLATION

Replacing the battery in the Transmitter:

The transmitter is operated by 1 x 3V CR2032 battery (battery included).

1. Remove the back plate of the transmitter by pushing in the tab at the bottom of the unit with a small flat screwdriver.
2. Remove the locking screw from the battery fastener of the battery compartment.
3. Replace the 3V Lithium (CR2032) battery with a battery of the same type and rating, ensuring the battery polarity matches the internal markings.
4. Replace the battery fastener with the locking screw to hold it in its compartment.
5. Replace the back plate of the transmitter.

Batteries are available on our site:
www.enquiries@matonguard.co.uk

INSTALLATION GUIDE

To avoid any undue interference from similar devices used by your neighbours, there is a unique selectable channel function button.

The transmission range of the wireless signal will be reduced if the signal is transmitted through a number of walls or ceilings.

Mounting the transmitter on metal surfaces will dramatically reduce the transmission range.

PORTABLE AND PLUG-IN CHIME RECEIVERS CAN BE COMBINED

PLUG-IN or PORTABLE CHIMES CAN BE USED WITH THE SAME PRESSURE MAT

DIGITAL CODE SETUP

The Care@Home-PMs is equipped with advanced random and self-learning operation codes.

First time operation:

After battery installation on the Transmitter and Receiver unit and plug the Plug-in receiver into any convenient wall socket, press the TEST button once to generate a new operation code. The chimes will sound to indicate the new operation code is registered and memorised.

If the Chimes sound with no reason, its most likely there is an interference from a nearby device:

1. Remove the battery from the Transmitter and Receiver and unplug the Plug-in receiver and rest for 10 seconds. This will erase the memory.
2. Insert the batteries and plug the plug-in receiver and the units are ready for the new code registration.

Important:

Repeat the above procedure after changing the Transmitter battery.

MELODY SELECTION

The melody selection is located inside the Transmitter unit. A small switch on the top right corner of the battery compartment marked with "S". Press the switch once, the melody will change. There are total 8 melodies for selection. The Plug-in Chime Receiver will play the Chime selected.

For any assistance or to purchase additional Chime Receiver units please contact:

Power Options Ltd. (Grand Offers) at:

enquiries@matonguard.co.uk.co.uk

Phone: 0208 660 6404

MAT-ON-GUARD

CARE & MOME-PMs



Model: PM



Your **Care@Home-PMs** product includes:

- ◆ A 72 x 39cm Electronic Pressure Mat
- ◆ A Transmitter device attached to the Electronic Pressure Mat
- ◆ Two Chime Receivers (Plug-in and Portable)
- ◆ 2 x AA Batteries for the Portable Receiver
- ◆ Portable and Plug-in receivers can be used combined with the same Mat
- ◆ Grip-Mat

INTRODUCTION

The MAT-ON-GUARD Care@Home-PMs product is simple to setup and easy to use.

Setting up:

1. Plug the Plug-in Receiver in an electric socket in the house.
2. Insert the batteries in the portable receiver.
3. Remove the Ribbon attached to the Transmitter to release the battery.
4. Press the TEST button on the Transmitter.
5. The Chime Receivers should chime.
6. Connect the Transmitter to the Electronic Pressure Mat.
7. Apply pressure on the Pressure Mat and the Chime Receivers should sound

The small radio transmitter which is connected to the pressure mat sends a signal to the receiver chime units indicating that the pressure mat is activated. The transmitter operating range is approximately 75 metres from the chime receivers, although that range may vary.

www.enquiries@matonguard.co.uk